

**SKIN CHRONICLE**  
Calibration Reference Card

10mm


Grid

0 10 20 30 40 50 60 70

skinchronicle.com

**HOW TO USE**

1. Place card flat next to skin lesion
2. Keep card and lesion at same distance from camera
3. Ensure ruler markings are clearly visible
4. Use grid or ruler edge for measurement calibration

**PRINTING INSTRUCTIONS**

Print at 100% scale (no fit-to-page) in landscape orientation.  
Cut along the corner guides. Card size: 85.6mm x 53.98mm (credit card size).